

Simulation report

Type of Freezer: Conventional/Radial freezer (e.g., CoolCell®)

Biomixture: DMSO + Culture medium

Mass fraction of DMSO: 0.01

Cooling rate (°C/min): 5




Figure 1: Cumulative mass distribution of ice growth velocity ($\mu\text{m/s}$)



Cryoconcentration (C/C_0)

Cryoconcentration is defined as the ratio of solute concentration after freezing (C) and initial solute concentration (C_0).




Figure 2: Cumulative mass distribution of cryoconcentration (C/C_0)

Freezing stress-time (min)

Freezing stress-time is defined as the time that a cell spends between the freezing temperature (T_f) and the vitrification temperature (T_g).




Figure 3: Cumulative mass distribution of freezing stress-time (min)



PRECISION CRYOSYSTEMS

***For more informations
or technical support:***

Please contact SMARTFREEZ by the
following e-mails:

info@smartfreez.com
support@smartfreez.com

Our engineers will get back to you.

www.smartfreez.com

